



# D-dimer Control Reference Value Table

## D-dimer Control

LOT 2025120801

 20251208

 20271207

Applicable Instruments	Unit	Concentration	Reference Value	Deviation Limit
CA1200 CA1203 CA1205 DC1210 DC1250 DC1280	μg/mL (FEU)	High value	1.79	±0.54

**【NOTE】**

1. The control should be stored in the refrigerator at 2 °C ~8 °C. After reconstitution, the control can be stable for 8 hours at room temperature (10 °C~30 °C), and can be stable for 12 hours at 2 °C~8 °C.
2. The control cannot be frozen after reconstitution.
3. The control is for IVD only.