


Table of Serum Amyloid A Protein Control (SAA) Reference Value

Serum Amyloid A Protein Control

LOT 2025121601

 2026-12-15

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DM79X、DM78X DM72X、DM71X (Applicable to technical file version A1.0 or higher)	Control 1	8.06	±2.42	DP-H10、DP-H12 DP-H13、DP-H15 DP-H16、DP-H17 (Applicable to technical file version A1.0 or higher)	Control 1	8.7	±2.6
	Control 2	63.39	±19.02		Control 2	53.5	±16.1