



# Coagulation Control Reference Value Table

## Coagulation Control-level 1

LOT 2025121901

 20251219

 20271218

Applicable Instruments	Parameter	Unit	Principle	Reference Value	Deviation Limit
CA1200 CA1203 CA1205 DC1210 DC1250 DC1280	PT (Lyophilized)	s	Optical method	12.5	±2.5
		s	Magnetic bead	12.5	±2.5
	PT (liquid)	s	Optical method	12.9	±2.6
		s	Magnetic bead	13.0	±2.6
	APTT (Lyophilized)	s	Optical method	30.0	±6.0
		s	Magnetic bead	29.3	±5.9
	APTT (liquid)	s	Optical method	33.1	±6.6
		s	Magnetic bead	32.0	±6.4
	TT (Lyophilized)	s	Optical method	17.4	±3.5
		s	Magnetic bead	18.8	±3.8
	TT (liquid)	s	Optical method	20.1	±4.0
		s	Magnetic bead	19.8	±4.0
	FIB (Lyophilized)	g/L	Optical method	2.86	±0.57
		g/L	Magnetic bead	2.94	±0.59
	FIB (liquid)	g/L	Optical method	2.73	±0.55
		g/L	Magnetic bead	2.75	±0.55
	AT-III (liquid)	%	Optical method	77.7	±15.5