


# Table of C-Reactive Protein Calibrator (CRP) Reference Value

## C-Reactive Protein Calibrator

LOT 2026010501

 2027-01-04

Applicable Instruments	DH51CRP, DH53CRP DH56CRP, D1-CRP D3-CRP, D5-CRP	DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP DH76CRP	DF1-CRP, DF3-CRP DF5-CRP, DF50CRP DF52CRP, DF53CRP
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	5.3	4.6	/	4.6	4.3
Calibrator 2	19.2	18.4	/	18.4	17.8
Calibrator 3	50.6	52.0	/	49.2	47.3
Calibrator 4	85.3	81.1	/	87.3	90.7
Calibrator 5	106.4	95.9	/	103.5	131.9
Calibrator 6	163.3	178.8	/	148.1	167.4
Applicable Instruments	DP-H10, DP-H12 DP-H13, DP-H15 DP-H16, DP-H17	DM79X, DM78X DM72X, DM71X	DH36X DH33X DH31X	DH-800 CRP series DH-800 CS series	DH-2000T series DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	5.2	5.32	4.5	3.94	4.42
Calibrator 2	22.1	19.38	22.1	17.94	16.93
Calibrator 3	56.8	50.65	49.0	48.25	47.19
Calibrator 4	96.6	85.04	72.8	79.19	71.78
Calibrator 5	114.0	103.83	95.8	91.31	81.26
Calibrator 6	190.0	165.61	153.8	156.57	151.85