


Table of C-Reactive Protein Control (CRP) Reference Value

C-Reactive Protein Control

LOT 2026010501

 2027-01-04

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP DH56CRP, D1-CRP D3-CRP, D5-CRP	Control 1	3.3	±1.0	DP-H10, DP-H12 DP-H13, DP-H15 DP-H16, DP-H17	Control 1	3.6	±1.1
	Control 2	14.3	±4.3		Control 2	15.8	±4.7
	Control 3	38.1	±11.4		Control 3	48.0	±14.4
DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	Control 1	2.9	±0.87	DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	Control 1	/	/
	Control 2	13.8	±4.13		Control 2	/	/
	Control 3	39.2	±11.75		Control 3	/	/
D6-CRP DH76CRP	Control 1	2.9	±0.9	DF1-CRP, DF3-CRP DF5-CRP, DF50CRP DF52CRP, DF53CRP	Control 1	2.7	±0.8
	Control 2	13.2	±4.0		Control 2	12.6	±3.8
	Control 3	37.2	±11.2		Control 3	38.6	±11.6
DM79X, DM78X DM72X, DM71X	Control 1	3.38	±1.01	DH36X DH33X DH31X	Control 1	2.8	±0.8
	Control 2	14.38	±4.31		Control 2	13.6	±4.1
	Control 3	38.78	±11.63		Control 3	39.5	±11.9
DH-2000T series DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T	Control 1	2.91	±0.87	DH-800 CRP series DH-800 CS series	Control 1	2.50	±0.75
	Control 2	12.78	±3.83		Control 2	13.06	±3.92
	Control 3	36.47	±10.94		Control 3	41.65	±12.50