

# D-dimer Calibrator Reference Value Table

## D-dimer Calibrator

[LOT] 2025122401

[M] 20251224

[E] 20271223

Applicable Instruments	Unit	Reference Value
CA1200 CA1203 CA1205 DC1210 DC1250 DC1280	µg/mL (FEU)	3.88

**【NOTE】**

1. The calibrator should be stored in the refrigerator at 2 °C~8 °C. After reconstitution, the calibrator can be stable for 8 hours at room temperature (10 °C~30 °C), and can be stable for 12 hours at 2 °C~8 °C.
2. The calibrator cannot be frozen after reconstitution.
3. The calibrator is for IVD only.