

	文件类型	技术文件	文件版本	V12.3.1	
	项目代号	/	页码	第0页 共4页	
	文件编号	GL DM-56-02-en			
文件名称： <p style="text-align: center;">CALIBRATOR NK0326 2026-04-10</p> 适用范围： 公司所有五分类仪器					
编制：石丽岐 日期：2026.03.03		审核：冯会娟 日期：2026.03.03		批准：冯会娟 日期：2026.03.03	
修 订 记 录					
版本	TCN/ECR/PCN编号	修订内容概要	修改人	批准人	批准日期
V12.3.0	/	新建	刘金盈	冯会娟	2026.02.12
V12.3.1	/	新增7003项目参考值	石丽岐	冯会娟	2026.03.03
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检				
由工程					

DM-CAL PLUS

HEMATOLOGY CALIBRATOR

Reference Values provided by DYMIND
CALIBRATOR
LOT NK0326

 2026-02-10

 2026-04-10

Applicable Instruments	Parameter	Unit	Reference Value
D1-CRP, D3-CRP, D5-CRP DH51CRP, DH53CRP, DH56CRP (Technical File Version A11.11 or higher)	WBC	$\times 10^9/L$	9.68
	RBC	$\times 10^{12}/L$	4.33
	HGB	g/L	131
	MCV	fL	89.7
	PLT	$\times 10^9/L$	221
UN71, UN73, UN76, DH73 (Technical File Version A6.3 or higher)	WBC	$\times 10^9/L$	9.36
	RBC	$\times 10^{12}/L$	4.43
	HGB	g/L	134
	MCV	fL	89.7
	PLT	$\times 10^9/L$	209
DH73 Vet (Technical File Version B5.5 or higher)	RBC	$\times 10^{12}/L$	4.50
	HGB	g/L	132
	MCV	fL	92.1
	PLT	$\times 10^9/L$	210
DF50, DF51, DF52 DF53, DF55, DF56 (Technical File Version A12.0 or higher)	WBC	$\times 10^9/L$	9.06
	RBC	$\times 10^{12}/L$	4.45
	HGB	g/L	131
	MCV	fL	87.5
	PLT	$\times 10^9/L$	224


【NOTE】


- The calibrator should be stored in refrigerator (2℃~8℃). After opening, it will keep stable for 7 days when it is stored airtight at 2℃~8℃.
- Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15℃~30℃).
- Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
- After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.

DM-CAL PLUS

HEMATOLOGY CALIBRATOR

Reference Values provided by DYMIND
CALIBRATOR
 NK0326

 2026-02-10

 2026-04-10



Applicable Instruments	Parameter	Unit	Reference Value
DF50, DF51, DF52 DF53, DF55, DF56 (Technical File Version B1.0 or higher)	WBC	$\times 10^9/L$	9.12
	RBC	$\times 10^{12}/L$	4.32
	HGB	g/L	132
	MCV	fL	88.1
	PLT	$\times 10^9/L$	209
DH71, DH76 (Technical File Version A4.1 to A10.3) DH51, DH53, DH56 (Technical File Version A8.1 or higher)	WBC	$\times 10^9/L$	9.20
	RBC	$\times 10^{12}/L$	4.26
	HGB	g/L	132
	MCV	fL	92.3
	PLT	$\times 10^9/L$	217
DH71, DH76 (Technical File Version A10.5 or higher)	WBC	$\times 10^9/L$	9.55
	RBC	$\times 10^{12}/L$	4.48
	HGB	g/L	131
	MCV	fL	87.7
	PLT	$\times 10^9/L$	221
D2-CRP, D7-CRP DH71CRP, DH73CRP (Technical File Version A6.0 to A6.5)	WBC	$\times 10^9/L$	9.42
	RBC	$\times 10^{12}/L$	4.42
	HGB	g/L	133
	MCV	fL	88.5
	PLT	$\times 10^9/L$	219

【NOTE】

- The calibrator should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 7 days when it is stored airtight at 2°C~8°C.
- Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
- Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
- After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.

DM-CAL PLUS

HEMATOLOGY CALIBRATOR

Reference Values provided by DYMIND
CALIBRATOR
 **NK0326**
 **2026-02-10**
 **2026-04-10**

Applicable Instruments	Parameter	Unit	Reference Value
D2-CRP, D7-CRP DH71CRP, DH73CRP (Technical File Version A6.6 or higher)	WBC	$\times 10^9/L$	9.69
	RBC	$\times 10^{12}/L$	4.37
	HGB	g/L	134
	MCV	fL	88.1
	PLT	$\times 10^9/L$	215
D6-CRP, DH76CRP (Technical File Version B2.2 or higher)	WBC	$\times 10^9/L$	9.60
	RBC	$\times 10^{12}/L$	4.52
	HGB	g/L	134
	MCV	fL	89.0
	PLT	$\times 10^9/L$	222
DF50 Vet, DF52 Vet DF55 Vet, DF56 Vet (Technical File Version A8.0 or higher and B1.0 or higher)	WBC	$\times 10^9/L$	8.58
	RBC	$\times 10^{12}/L$	4.36
	HGB	g/L	129
	MCV	fL	89.1
	PLT	$\times 10^9/L$	216


【NOTE】


1. The calibrator should be stored in refrigerator (2℃~8℃). After opening, it will keep stable for 7 days when it is stored airtight at 2℃~8℃.
2. Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15℃~30℃).
3. Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.

DM-CAL PLUS

HEMATOLOGY CALIBRATOR

Reference Values provided by DYMIND
CALIBRATOR
 NK0326

 2026-02-10

 2026-04-10

Applicable Instruments	Parameter	Unit	Reference Value
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Technical File Version A6.1 or higher)	WBC	$\times 10^9/L$	9.62
	RBC	$\times 10^{12}/L$	4.35
	HGB	g/L	131
	MCV	fL	88.9
	PLT	$\times 10^9/L$	222
DM71X, DM72X DM75X, DM77X DM78X, DM79X	WBC	$\times 10^9/L$	9.29
	RBC	$\times 10^{12}/L$	4.36
	HGB	g/L	132
	MCV	fL	87.9
	PLT	$\times 10^9/L$	216
UN71 Vet, UN73 Vet (Technical File Version A1.0 or higher)	WBC	$\times 10^9/L$	8.91
	RBC	$\times 10^{12}/L$	4.40
	HGB	g/L	133
	MCV	fL	89.0
	PLT	$\times 10^9/L$	219
DM60 Vet, DM61 Vet DM62 Vet, DM63 Vet (Technical File Version A1.4 or higher)	WBC	$\times 10^9/L$	7.67
	RBC	$\times 10^{12}/L$	3.91
	HGB	g/L	126
	MCV	fL	89.4
	PLT	$\times 10^9/L$	215

【NOTE】

- The calibrator should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 7 days when it is stored airtight at 2°C~8°C.
- Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
- Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
- After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.