


# Table of C-Reactive Protein Calibrator (CRP) Reference Value

## C-Reactive Protein Calibrator

LOT 2026031402

 2027-03-13

Applicable Instruments	DH51CRP, DH53CRP DH56CRP, D1-CRP D3-CRP, D5-CRP	DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP DH76CRP	DF1-CRP, DF3-CRP DF5-CRP, DF50CRP DF52CRP, DF53CRP
<b>Level (multi-point)</b>	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
<b>Calibrator 1</b>	5.1	4.5	4.9	4.5	4.9
<b>Calibrator 2</b>	19.4	20.7	20.2	21.8	19.7
<b>Calibrator 3</b>	43.8	46.0	46.0	44.1	44.2
<b>Calibrator 4</b>	82.8	88.4	87.0	74.4	85.8
<b>Calibrator 5</b>	104.5	104.3	103.4	99.1	102.8
<b>Calibrator 6</b>	177.2	161.2	155.8	146.0	148.3
Applicable Instruments	DP-H10, DP-H12 DP-H13, DP-H15 DP-H16, DP-H17	DM79X, DM78X DM72X, DM71X DM75X	DH36X DH33X DH31X	DH-800 CRP series DH-800 CS series	DH-2000T series DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T
<b>Level (multi-point)</b>	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
<b>Calibrator 1</b>	6.2	5.21	5.3	4.22	3.97
<b>Calibrator 2</b>	22.9	21.57	21.5	20.82	19.32
<b>Calibrator 3</b>	59.7	48.17	49.4	47.41	43.42
<b>Calibrator 4</b>	99.3	76.69	85.8	76.15	74.73
<b>Calibrator 5</b>	114.0	91.33	104.3	88.12	88.03
<b>Calibrator 6</b>	155.9	155.60	166.2	154.27	142.51

C-reactive Protein Calibrator (Lot 2026031402) are compatible with C-reactive Protein Reagent lots 2026031401 and 2026031402.