


# Table of C-Reactive Protein Control (CRP) Reference Value

## C-Reactive Protein Control

LOT 2026031402

 2027-03-13

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP DH56CRP, D1-CRP D3-CRP, D5-CRP	Control 1	<b>3.0</b>	±0.9	DP-H10, DP-H12 DP-H13, DP-H15 DP-H16, DP-H17	Control 1	<b>3.7</b>	±1.1
	Control 2	<b>13.3</b>	±4.0		Control 2	<b>16.7</b>	±5.0
	Control 3	<b>36.4</b>	±10.9		Control 3	<b>45.7</b>	±13.7
DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	Control 1	<b>2.9</b>	±0.9	DH71CRP, DH73CRP D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	Control 1	<b>3.1</b>	±0.9
	Control 2	<b>13.0</b>	±3.9		Control 2	<b>13.4</b>	±4.0
	Control 3	<b>36.9</b>	±11.1		Control 3	<b>36.3</b>	±10.9
D6-CRP DH76CRP	Control 1	<b>2.8</b>	±0.8	DF1-CRP, DF3-CRP DF5-CRP, DF50CRP DF52CRP, DF53CRP	Control 1	<b>2.9</b>	±0.9
	Control 2	<b>14.6</b>	±4.4		Control 2	<b>13.7</b>	±4.1
	Control 3	<b>37.9</b>	±11.4		Control 3	<b>36.6</b>	±11.0
DM79X, DM78X DM72X, DM71X DM75X	Control 1	<b>3.20</b>	±0.96	DH36X DH33X DH31X	Control 1	<b>3.4</b>	±1.0
	Control 2	<b>15.31</b>	±4.59		Control 2	<b>14.5</b>	±4.3
	Control 3	<b>40.12</b>	±12.04		Control 3	<b>38.6</b>	±11.6
DH-2000T series DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T	Control 1	<b>2.45</b>	±0.74	DH-800 CRP series DH-800 CS series	Control 1	<b>3.06</b>	±0.92
	Control 2	<b>12.95</b>	±3.89		Control 2	<b>15.62</b>	±4.69
	Control 3	<b>35.61</b>	±10.68		Control 3	<b>41.49</b>	±12.45

C-reactive Protein Control (Lot 2026031402) are compatible with C-reactive Protein Reagent lots 2026031401 and 2026031402.