



DM-BF

HEMATOLOGY CONTROL

Reference Values provided by DYMIND
CONTROL
 2026-02-10

 2026-05-10

Applicable Instruments	Parameter	Unit	LOT BF0326L	LOT BF0326N	LOT BF0326H
DYMIND DH-800 series DH-800 CRP series DH-800 CS series (Technical File Version A1.3 or higher)	WBC-BF	10 ⁹ /L	0.074 ±0.050	0.228 ±0.069	0.911 ±0.137
	TC-BF#	10 ⁹ /L	0.074 ±0.050	0.229 ±0.069	0.913 ±0.137
	MN#	10 ⁹ /L	0.025 ±0.015	0.070 ±0.014	0.146 ±0.088
	MN%	%	33.8 ±25.0	30.9 ±17.0	16.0 ±8.0
	PMN#	10 ⁹ /L	0.049 ±0.025	0.158 ±0.079	0.765 ±0.383
	PMN%	%	66.2 ±33.1	69.1 ±27.7	84.0 ±25.2
	RBC-BF	10 ¹² /L	0.024 ±0.020	0.071 ±0.060	0.480 ±0.100

【NOTE】

1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 30 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
- 3.Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



BF0326L



BF0326N



BF0326H